

CLIENT DETAILS

LABORATORY DETAILS

Contact **Brendan Stuart**
 Client **PREMISE**
 Address **LEVEL 1
 100 BRUNSWICK STREET
 FORTITUDE VALLEY
 NSW 2006**
 Telephone **61 2 6939 5000**
 Facsimile **(Not specified)**
 Email **Brendan.stuart@premise.com.au**
 Project **217500 - Lithgow SWF**
 Order Number **217500**
 Samples **1**

Manager **Shane McDermott**
 Laboratory **SGS Alexandria Environmental**
 Address **Unit 16, 33 Maddox St
 Alexandria NSW 2015**
 Telephone **+61 2 8594 0400**
 Facsimile **+61 2 8594 0499**
 Email **au.environmental.sydney@sgs.com**
 SGS Reference **SE287383 R0**
 Date Received **8/8/2025**
 Date Reported **15/8/2025**

COMMENTS

Accredited for compliance with ISO/IEC 17025 - Testing. NATA accredited laboratory 2562(4354).

SIGNATORIES

Gemma MOON
Organic Chemist

Kamrul AHSAN
Senior Chemist

Ly Kim HA
Organic Section Head

Shane MCDERMOTT
Laboratory Manager

Ying Ying ZHANG
Laboratory Technician

Volatile Petroleum Hydrocarbons in Water [AN433] Tested: 8/8/2025

			SW1
			WATER
			-
			7/8/2025
PARAMETER	UOM	LOR	SE287383.001
Benzene (F0)	µg/L	0.5	<0.5
TRH C6-C9	µg/L	40	<40
TRH C6-C10	µg/L	50	<50
TRH C6-C10 minus BTEX (F1)	µg/L	50	<50

TRH (Total Recoverable Hydrocarbons) in Water [AN403] Tested: 11/8/2025

			SW1
			WATER
			-
			7/8/2025
PARAMETER	UOM	LOR	SE287383.001
TRH C10-C14	µg/L	50	<50
TRH C15-C28	µg/L	200	220
TRH C29-C36	µg/L	200	<200
TRH C37-C40	µg/L	200	<200
TRH >C10-C16	µg/L	60	<60
TRH >C10-C16 - Naphthalene (F2)	µg/L	60	<60
TRH >C16-C34 (F3)	µg/L	500	<500
TRH >C34-C40 (F4)	µg/L	500	<500
TRH C10-C40	µg/L	320	<320

OC Pesticides in Water [AN420] Tested: 11/8/2025

PARAMETER	UOM	LOR	SW1
			WATER - 7/8/2025 SE287383.001
Hexachlorobenzene (HCB)	µg/L	0.1	<0.1
Alpha BHC	µg/L	0.1	<0.1
Lindane (gamma BHC)	µg/L	0.1	<0.1
Heptachlor	µg/L	0.1	<0.1
Aldrin	µg/L	0.1	<0.1
Beta BHC	µg/L	0.1	<0.1
Delta BHC	µg/L	0.1	<0.1
Heptachlor epoxide	µg/L	0.1	<0.1
Alpha Endosulfan	µg/L	0.1	<0.1
Gamma Chlordane	µg/L	0.1	<0.1
Alpha Chlordane	µg/L	0.1	<0.1
p,p'-DDE	µg/L	0.1	<0.1
Dieldrin	µg/L	0.1	<0.1
Endrin	µg/L	0.1	<0.1
Beta Endosulfan	µg/L	0.1	<0.1
p,p'-DDD	µg/L	0.1	<0.1
p,p'-DDT	µg/L	0.1	<0.1
Endosulfan sulphate	µg/L	0.1	<0.1
Endrin aldehyde	µg/L	0.1	<0.1
Methoxychlor	µg/L	0.1	<0.1
Endrin ketone	µg/L	0.1	<0.1
Mirex	µg/L	0.1	<0.1
Total OC	µg/L	1	<1
Total OC	µg/L	1	<1

OP Pesticides in Water [AN420] Tested: 11/8/2025

PARAMETER	UOM	LOR	SW1
			WATER - 7/8/2025 SE287383.001
Dichlorvos	µg/L	0.5	<0.5
Dimethoate	µg/L	0.5	<0.5
Diazinon (Dimpylate)	µg/L	0.5	<0.5
Fenitrothion	µg/L	0.2	<0.2
Malathion	µg/L	0.2	<0.2
Chlorpyrifos (Chlorpyrifos Ethyl)	µg/L	0.2	<0.2
Parathion-ethyl (Parathion)	µg/L	0.2	<0.2
Bromophos Ethyl	µg/L	0.2	<0.2
Methidathion	µg/L	0.5	<0.5
Ethion	µg/L	0.2	<0.2
Azinphos-methyl	µg/L	0.2	<0.2

Total Phenolics in Water [AN295] Tested: 13/8/2025

			SW1
			WATER
			-
			7/8/2025
PARAMETER	UOM	LOR	SE287383.001
Total Phenols	mg/L	0.05	<0.05

pH in water [AN101] Tested: 8/8/2025

			SW1
			WATER
			-
			7/8/2025
PARAMETER	UOM	LOR	SE287383.001
pH**	No unit	-	7.4

Conductivity and TDS by Calculation - Water [AN106] Tested: 8/8/2025

			SW1
			WATER
			-
			7/8/2025
PARAMETER	UOM	LOR	SE287383.001
Conductivity @ 25 C	µS/cm	2	500

Ammonia Nitrogen by Discrete Analyser [AN291] Tested: 8/8/2025

			SW1
			WATER
			-
			7/8/2025
PARAMETER	UOM	LOR	SE287383.001
Ammonia Nitrogen, NH ₃ as N	mg/L	0.01	8.9

Anions by Ion Chromatography in Water [AN245] Tested: 11/8/2025

			SW1
			WATER
			-
			7/8/2025
PARAMETER	UOM	LOR	SE287383.001
Fluoride	mg/L	0.1	0.12
Chloride	mg/L	1	37
Sulfate, SO4	mg/L	1	38
Nitrate Nitrogen, NO3-N	mg/L	0.005	3.8

Alkalinity [AN135] Tested: 12/8/2025

			SW1
			WATER
			-
			7/8/2025
PARAMETER	UOM	LOR	SE287383.001
Total Alkalinity as CaCO3	mg/L	5	140

Total Phosphorus by Kjeldahl Digestion DA in Water [AN279/AN293] Tested: 11/8/2025

			SW1
			WATER
			-
			7/8/2025
PARAMETER	UOM	LOR	SE287383.001
Total Phosphorus (Kjeldahl Digestion) as P	mg/L	0.02	0.39

Total and Volatile Suspended Solids (TSS / VSS) [AN114] Tested: 11/8/2025

			SW1
			WATER
			-
			7/8/2025
PARAMETER	UOM	LOR	SE287383.001
Total Suspended Solids Dried at 103-105°C	mg/L	5	220

COD in Water [AN179/AN181] Tested: 13/8/2025

			SW1
			WATER
			-
			7/8/2025
PARAMETER	UOM	LOR	SE287383.001
Chemical Oxygen Demand	mg/L	10	80

Forms of Carbon [AN190] Tested: 11/8/2025

			SW1
			WATER
			-
			7/8/2025
PARAMETER	UOM	LOR	SE287383.001
Total Organic Carbon as NPOC	mg/L	0.2	18

Trace Metals (Dissolved) in Water by ICPMS [AN318] Tested: 13/8/2025

			SW1
			WATER
			-
			7/8/2025
			SE287383.001
PARAMETER	UOM	LOR	
Iron	µg/L	5	840
Manganese	µg/L	1	5
Aluminium	µg/L	5	870

Metals in Water (Dissolved) by ICPOES [AN320] Tested: 14/8/2025

			SW1
			WATER
			-
			7/8/2025
			SE287383.001
PARAMETER	UOM	LOR	
Calcium, Ca	mg/L	0.2	34
Magnesium, Mg	mg/L	0.1	10
Sodium, Na	mg/L	0.5	29
Potassium, K	mg/L	0.1	22

Trace Metals (Total) in Water by ICPMS [AN022/AN318] Tested: 13/8/2025

			SW1
			WATER
			-
			7/8/2025
PARAMETER	UOM	LOR	SE287383.001
Total Chromium	µg/L	1	6

Hexavalent Chromium in water by Discrete Analyser [AN283] Tested: 8/8/2025

			SW1
			WATER
			-
			7/8/2025
PARAMETER	UOM	LOR	SE287383.001
Hexavalent Chromium, Cr6+	mg/L	0.004	<0.004

METHOD

METHODOLOGY SUMMARY

- AN020** Unpreserved water sample is filtered through a 0.45µm membrane filter and acidified with nitric acid similar to APHA3030B.
- AN022/AN318** Following acid digestion of un filtered sample, determination of elements at trace level in waters by ICP-MS technique, referenced to USEPA 6020B and USEPA 200.8 (5.4).
- AN022** The water sample is digested with Nitric Acid and made up to the original volume similar to APHA3030E.
- AN101** pH in Soil Sludge Sediment and Water: pH is measured electrometrically using a combination electrode (glass plus reference electrode) and is calibrated against 3 buffers purchased commercially. For soils, an extract with water is made at a ratio of 1:5 and the pH determined and reported on the extract. Reference APHA 4500-H+.
- AN106** Conductivity and TDS by Calculation: Conductivity is measured by meter with temperature compensation and is calibrated against a standard solution of potassium chloride. Conductivity is generally reported as µmhos/cm or µS/cm @ 25°C. For soils, an extract with water is made at a ratio of 1:5 and the EC determined and reported on the extract, or calculated back to the as-received sample. Total Dissolved Salts can be estimated from conductivity using a conversion factor, which for natural waters, is in the range 0.55 to 0.75. SGS use 0.6. Reference APHA 2510 B.
- AN106** Salinity may be calculated in terms of NaCl from the sample conductivity. This assumes all soluble salts present, measured by the conductivity, are present as NaCl.
- AN114** Total Suspended and Volatile Suspended Solids: The sample is homogenised by shaking and a known volume is filtered through a pre-weighed GF/C filter paper and washed well with deionised water. The filter paper is dried and reweighed. The TSS is the residue retained by the filter per unit volume of sample. Reference APHA 2540 D. Internal Reference AN114
- AN135** Alkalinity (and forms of) by Titration: The sample is titrated with standard acid to pH 8.3 (P titre) and pH 4.5 (T titre) and permanent and/or total alkalinity calculated. The results are expressed as equivalents of calcium carbonate or recalculated as bicarbonate, carbonate and hydroxide. Reference APHA 2320. Internal Reference AN135
- AN181** Analysis of COD by Semi Closed Reflux: The sample is refluxed with strong acid and a known excess of oxidant. After digestion the unreduced oxidant is back titrated to determine the amount of oxidant consumed. The chemically oxidised matter is calculated in terms of oxygen equivalents. Reference APHA 5220 B.
- AN190** TOC and DOC in Water: A homogenised micro portion of sample is injected into a heated reaction chamber packed with an oxidative catalyst that converts organic carbon to carbon dioxide. The CO₂ is measured using a non-dispersive infrared detector. The process is fully automated in a commercially available analyser. If required a sugar value can be calculated from the TOC result. Reference APHA 5310 B.
- AN190** Chemical oxygen demand can be calculated/estimated based on the O₂/C relation as 2.67*NPOC (TOC). This is an estimate only and the factor will vary with sample matrix so results should be interpreted with caution.
- AN245** Anions by Ion Chromatography: A water sample is injected into an eluent stream that passes through the ion chromatographic system where the anions of interest ie Br, Cl, NO₂, NO₃ and SO₄ are separated on their relative affinities for the active sites on the column packing material. Changes to the conductivity and the UV-visible absorbance of the eluent enable identification and quantitation of the anions based on their retention time and peak height or area. APHA 4110 B
- AN279/AN293(Sydney)** The sample is digested with Sulphuric acid, K₂SO₄ and CuSO₄. All forms of phosphorus are converted into orthophosphate. The digest is cooled and placed on the discrete analyser for colorimetric analysis.
- AN283** Hexavalent Chromium via DA: Soluble hexavalent chromium forms a red/violet colour with diphenylcarbazide in acidic solution. This procedure is very sensitive and nearly specific for Cr⁶⁺. If total chromium is also measured the trivalent form of chromium Cr³⁺ can be calculated from the difference (Total Cr - Cr⁶⁺). Reference APHA3500CrB.
- AN291** Ammonia in solution reacts with hypochlorite ions from Sodium Dichloroisocyanate, and salicylate in the presence of Sodium Nitroprusside to form indophenol blue and measured at 660 nm by Discrete Analyser.
- AN295** The water sample or extract of sample is distilled in a phosphoric acid stream. Phenolic compounds in the distillate react with a reagent stream of potassium hexacyanoferrate(III) and 4-Amino-2,3-dimethyl-3-pyrazolin-5-one in an alkaline medium to form a coloured complex which is analysed spectrophotometrically onboard a continuous flow analyser.
- AN318** Determination of elements at trace level in waters by ICP-MS technique,, referenced to USEPA 6020B and USEPA 200.8 (5.4).

<p>AN320</p>	<p>Metals by ICP-OES: Samples are preserved with 10% nitric acid for a wide range of metals and some non-metals. This solution is measured by Inductively Coupled Plasma. Solutions are aspirated into an argon plasma at 8000-10000K and emit characteristic energy or light as a result of electron transitions through unique energy levels. The emitted light is focused onto a diffraction grating where it is separated into components .</p>
<p>AN320</p>	<p>Photomultipliers or CCDs are used to measure the light intensity at specific wavelengths. This intensity is directly proportional to concentration. Corrections are required to compensate for spectral overlap between elements . Reference APHA 3120 B.</p>
<p>AN403</p>	<p>Total Recoverable Hydrocarbons: Determination of Hydrocarbons by gas chromatography after a solvent extraction. Detection is by flame ionisation detector (FID) that produces an electronic signal in proportion to the combustible matter passing through it. Total Recoverable Hydrocarbons (TRH) are routinely reported as four alkane groupings based on the carbon chain length of the compounds: C6-C9, C10-C14, C15-C28 and C29-C36 and in recognition of the NEPM 1999 (2013), >C10-C16 (F2), >C16-C34 (F3) and >C34-C40 (F4). Where F2 is corrected for Naphthalene, the VOC data for Naphthalene is used.</p>
<p>AN403</p>	<p>Additionally, the volatile C6-C9/C6-C10 fractions may be determined by a purge and trap technique and GC/MS because of the potential for volatiles loss. Total Recoverable Hydrocarbons - Silica (TRH-Silica) follows the same method of analysis after silica gel cleanup of the solvent extract. Aliphatic/Aromatic Speciation follows the same method of analysis after fractionation of the solvent extract over silica with differential polarity of the eluent solvents.</p>
<p>AN403</p>	<p>The GC/FID method is not well suited to the analysis of refined high boiling point materials (ie lubricating oils or greases) but is particularly suited for measuring diesel, kerosene and petrol if care to control volatility is taken . This method will detect naturally occurring hydrocarbons, lipids, animal fats, phenols and PAHs if they are present at sufficient levels, dependent on the use of specific cleanup/fractionation techniques. Reference USEPA 3510B, 8015B.</p>
<p>AN420</p>	<p>SVOC Compounds: Semi-Volatile Organic Compounds (SVOCs) including OC, OP, PCB, Herbicides, PAH, Phthalates and Speciated Phenols in soils, sediments and waters are determined by GCMS/ECD technique following appropriate solvent extraction process (Based on USEPA 3500C and 8270D).</p>
<p>AN433</p>	<p>VOCs and C6-C9 Hydrocarbons by GC-MS P&T: VOC's are volatile organic compounds. The sample is presented to a gas chromatograph via a purge and trap (P&T) concentrator and autosampler and is detected with a Mass Spectrometer (MSD). Solid samples are initially extracted with methanol whilst liquid samples are processed directly. References: USEPA 5030B, 8020A, 8260.</p>
<p>Calculation</p>	<p>Free and Total Carbon Dioxide may be calculated using alkalinity forms only when the samples TDS is <500mg/L. If TDS is >500mg/L free or total carbon dioxide cannot be reported . APHA4500CO2 D.</p>

FOOTNOTES

*	NATA accreditation does not cover the performance of this service.	-	Not analysed.	UOM	Unit of Measure.
**	Indicative data, theoretical holding time exceeded.	NVL	Not validated.	LOR	Limit of Reporting.
***	Indicates that both * and ** apply.	IS	Insufficient sample for analysis.	↑↓	Raised/lowered Limit of Reporting.
		LNR	Sample listed, but not received.		

Unless it is reported that sampling has been performed by SGS, the samples have been analysed as received. Solid samples expressed on a dry weight basis.

Where "Total" analyte groups are reported (for example, Total PAHs, Total OC Pesticides) the total will be calculated as the sum of the individual analytes, with those analytes that are reported as <LOR being assumed to be zero. The summed (Total) limit of reporting is calculated by summing the individual analyte LORs and dividing by two. For example, where 16 individual analytes are being summed and each has an LOR of 0.1 mg/kg, the "Totals" LOR will be 1.6 / 2 (0.8 mg/kg). Where only 2 analytes are being summed, the "Total" LOR will be the sum of those two LORs.

Some totals may not appear to add up because the total is rounded after adding up the raw values.

If reported, measurement uncertainty follow the ± sign after the analytical result and is expressed as the expanded uncertainty calculated using a coverage factor of 2, providing a level of confidence of approximately 95%, unless stated otherwise in the comments section of this report.

Results reported for samples tested under test methods with codes starting with ARS-SOP, radionuclide or gross radioactivity concentrations are expressed in becquerel (Bq) per unit of mass or volume or per wipe as stated on the report. Becquerel is the SI unit for activity and equals one nuclear transformation per second.

Note that in terms of units of radioactivity:

- a. 1 Bq is equivalent to 27 pCi
- b. 37 MBq is equivalent to 1 mCi

For results reported for samples tested under test methods with codes starting with ARS-SOP, less than (<) values indicate the detection limit for each radionuclide or parameter for the measurement system used. The respective detection limits have been calculated in accordance with ISO 11929.

The QC and MU criteria are subject to internal review according to the SGS QAQC plan and may be provided on request or alternatively can be found here: www.sgs.com.au/en-gb/environment-health-and-safety.

This document is issued by the Company under its General Conditions of Service accessible at www.sgs.com/en/Terms-and-Conditions.aspx. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein.

Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client only. Any unauthorized alteration, forgery or